TRACOE softpad Tracheostomy Stoma Pad

TRACOE softpad Tracheostomapolster

REF 969

EN	Instructions for Use
DE	Gebrauchsanweisung
NL	Gebruiksaanwijzing
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
SV	Bruksanvisning
DA	Brugsanvisning
NO	Bruksanvisning
FI	Käyttöohje
PL	Instrukcja używania

EN / Warning · **DE** / Achtung · **NL** / Let op · **FR** / Attention · **IT** / Attenzione · **ES** / Atención · **PT** / Atenção · **SV** / OBS · **DA** / Vigtigt · **NO** / OBS · **FI** / Huomio · **PL** / Ostrzeżenie ·

EN / Consult Instructions for Use · **DE** / Gebrauchsanweisung beachten · **NL** / Rekening houden met de gebruiksaanwijzing · **FR** / Respecter le mode d'emploi · **IT** / Consultare le istruzioni per l'uso · **ES** / Siga las instrucciones de uso · **PT** / Ter atenção às instruções de utilização · **SV** / Se bruksanvisningen · **DA** / Følg brugsanvisningen · **NO** / Følg bruksanvisningen · **FI** / Noudata käyttöohjetta · **PL** / Zazrzyj do instrukcji używania ·

EN / Medical Device · **DE** / Medizinprodukt · **NL** / Medisch hulpmiddel · **FR** / Dispositif médical · **IT** / Dispositivo medico · **ES** / Dispositivo médico · **PT** / Dispositivo médico · **SV** / Medicinteknisk produkt · **DA** / Medicinsk udstyr · **NO** / Medisinsk enhet · **FI** / Lääkinnällinen laite · **PL** / Wyrób medyczny ·

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TRACOE softpad Tracheostomy Stoma Pad

REF 969

Note: Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your patients' and your own safety, please observe the following safety information.

1. Intended Use and Indications for Use

The TRACOE softpad tracheostomy stoma pad is used to provide padding between the tracheostomy tube and the tracheostoma. It is intended for single use only.

Clinical Benefit: The TRACOE softpad tracheostomy stoma pad is used to provide padding between the tracheostomy tube and the tracheostoma. It helps to absorb secretions and protect the skin beneath the neck flange.

Patient Population: The product is intended for children, adolescents (≥ 12 - 21 years), and adults.

Clinical Use: The product is intended for mechanically ventilated and self-breathing patients in hospitals, pre-hospitals (EMS), extended care facilities, or outpatient clinics, or home care.

Intended User: The product can be used by medical staff trained in tracheostomy care or individuals trained by professionals.

Indications for Use: Use of the TRACOE softpad tracheostomy stoma pad is indicated in patients whose tracheostomas have healed and are epithelialized.

Single Use and Useful Life: The TRACOE softpad tracheostomy stoma pad is single use only.

The TRACOE softpad tracheostomy stoma pad should be replaced at least on a daily basis. If there is increased secretion from the tracheostoma, more frequent replacement may be necessary.

Caution:

A prolonged and repeated use of a single tracheostomy stoma pad may result in material safety and biocompatibility issues.

2. General Description

The TRACOE softpad tracheostomy stoma pad is made of absorbent open-pore polyurethane foam. With its circular opening and notches, the product fits outer cannula sizes 5 to 10.

The tracheostomy stoma pad is oval, 58 mm wide and 36 mm high with a thickness of 5 mm; in the center it has a circular opening of 10 mm in diameter with four notches arranged around the sides at 90° to each other.

The TRACOE softpad tracheostomy stoma pad can be used in combination with other compresses for tracheostomy care.

3. Contraindications

- The tracheostomy stoma pad must not be used for sterile wound treatment.
- The tracheostomy stoma pad must not be applied on fresh and unhealed wounds.
- If secretion is heavy, the tracheostomy stoma pad should not be used. In this case, larger and/or highly absorbent compresses (e.g. TRACOE purofoam Tracheal Compress REF 958/REF 959) should be used.

4. General Precautions

- Follow the instructions for use of the respective tracheostomy tube in use. Contact the manufacturer if there are any questions, or if assistance is required.
- The product should be inspected for integrity and function prior to use. Verify that the material is not brittle or torn, and that there is an absence of kinks, tears or cuts. If the product is damaged, it should be replaced with a new product.
- Improper storage conditions may result in product damage.
- Do not perform a tracheostomy tube change in favor of exchanging the TRACOE softpad tracheostomy stoma pad. In this case, tearing and replacement with the tube in situ is preferable.

5. Warnings

- Reprocessing (including re-sterilisation) is not allowed, this may influence the material and function of the product. The products are single use only.
- Modifications of TRACOE products - except for the modifications explicitly described in the instructions for use (slitting) - are not allowed. TRACOE will not be responsible for modified products.

- In the event of allergic skin reactions such as itching, redness, rash or other signs of intolerance, the TRACOE softpad tracheostomy stoma pad must be removed, and a medical professional should be contacted immediately.
- When using sterile tracheostomy tubes, it is recommended that the TRACOE softpad tracheostomy stoma pad is already slit by tearing before its application, so that the application can take place without contact with the distal end of the tracheostomy tube. Sterility of the tracheostomy tube is no longer guaranteed once the stoma pad is pulled over the distal end of the tube.
- When applying or removing the TRACOE softpad tracheostomy stoma pad, make sure that the tracheostomy tube does not slip out of place.

6. Side Effects

Typical side effects of the use of stoma pads include infection, irritation and/or damage of the surrounding skin, the tracheostoma and the airway, pressure marks, discomfort, and local or systemic allergic reaction to the pad material. In case of an adverse event please contact a medical professional immediately.

7. Functional Description

1. Inspect the packaging to ensure it is secure and undamaged.
2. Open the package and inspect the device for damages prior to use.
3. Verify that the material is not brittle or torn and there are no kinks, tears or cuts.


Image 1

4. Before applying the TRACOE softpad tracheostomy stoma pad, a slit can be adjusted to the patient's needs. Slitting of the stoma pad should be performed by tearing at one of the pre-slit points from the inside to the outside, **see Image 1**. Afterwards, it should be ensured that there are no loose particles on the stoma pad, as these could get into the stoma or be aspirated.


Image 2

5. The tracheostomy stoma pad can be placed with each side on the patient's skin. It should be positioned between the tracheostomy tube and the skin, **see Image 2**. When using the tracheostomy stoma pad without slitting it, the tracheostomy tube is inserted through the opening in the center of the pad.

Caution:

Sterility of the tracheostomy tube is no longer guaranteed once the stoma pad is pulled over the distal end of the tube.

6. Option 1 (not pre-slitted): The tracheostomy tube can be inserted into the stoma according to the instructions for use of the manufacturer with the TRACOE softpad tracheostomy stoma pad in place.

Option 2 (pre-slitted): The tracheostomy tube is already inserted into the stoma according to the instructions for use of the manufacturer. Then, the TRACOE softpad tracheostomy stoma pad is gently slid between the neck plate and the patient's skin.

Caution:

When applying or removing the TRACOE softpad tracheostomy stoma pad, make sure that the tracheostomy tube does not slip out of place.

7. The TRACOE softpad tracheostomy stoma pad should be replaced at least on a daily basis. If there is increased secretion from the tracheostoma, more frequent replacement may be necessary. If the stoma pad needs to be changed at least once a day or several times a day, a slit should be applied.

Remove the TRACOE softpad tracheostomy stoma pad either when withdrawing the tracheostomy tube for cleaning or by tearing the stoma pad by hand, if it has not been slit prior to tube insertion.

Caution:

Do not perform a tracheostomy tube change in favor of exchanging the TRACOE softpad tracheostomy stoma pad. In this case, tearing and replacement with the tube in situ is preferable.

Caution:

Make sure that there are no loose particles left around the tracheostoma when removing the TRACOE softpad tracheostomy stoma pad.

8. Storage

Store the TRACOE products in their original packaging according to the conditions displayed on the packaging.

9. Packaging

The product is provided non-sterile and is individually packed. TRACOE softpad tracheostomy stoma pads do not require a sterile environment during normal use.

10. Disposal

Used products are to be disposed of in accordance with national regulations, waste management plans, or clinical procedures governing biohazardous waste materials, e.g. the direct disposal in a tear and moisture-resistant and secure bag or container, which is routed to the local waste disposal system for contaminated medical products.

For further recommendations, contact your hygiene officer in health facilities, or the local waste management for homecare use.

11. Returns and Complaints

Returned products, that have been used, will only be accepted, if TRACOE has agreed to the return and a completed decontamination certificate with complaint report is enclosed with the device. These forms are available either directly from TRACOE medical, or via the website www.tracoe.com.

If the device is involved in a reportable incident, as defined in local medical device legislation, please contact TRACOE medical (complaints@tracoe.com), and the appropriate regulatory body in the country of use.

12. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be affected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either from TRACOE medical GmbH or on our website at www.tracoe.com.

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