



TheraBite®

ActiveBand™



1)


2 b)


2 c)


3 a)


3 b)


3 c)


3 d)


3 e)


3 f)

4 a)

4 b)

4 c)

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Contents

EN - ENGLISH	7
DE - DEUTSCH	10
NL - NEDERLANDS.....	13
FR - FRANÇAIS.....	17
IT - ITALIANO.....	21
ES - ESPAÑOL.....	25
PT - PORTUGUÊS	28
SV - SVENSKA.....	32
DA - DANSK	35
NO - NORSK.....	39
FI -SUOMI	42
PL - POLSKI	45
RU - РУССКИЙ	49
JA - 日本語	53
ZHTW - 繁體中文.....	56
HE - العربية.....	64
AR - עברית	68

ENGLISH

Before starting any rehabilitation activity involving the TheraBite system and ActiveBand every user should consider the benefits of medically supervised programs. Clinicians may select and combine exercises to form programs that are appropriate for each user's abilities and goals.

1. Descriptive information

1.1 Intended Use

The TheraBite ActiveBand is an elastic silicone band used together with the TheraBite Jaw Mobilizer to increase and/or maintain muscle strength and endurance of the muscles of mastication (chewing muscles).

The TheraBite ActiveBand is intended for single-patient use only.

1.2 Contraindications

Contraindications as described in the instructions for use for TheraBite Jaw Motion Rehabilitation System apply.

1.3 Description of the Device

TheraBite ActiveBand is an unsterile elastic silicone band for single patient use. It is looped around the Lever of the TheraBite Jaw Mobilizer (included in the product TheraBite Jaw Motion Rehabilitation System). The ActiveBand's position on the Load Scale determines the level of resistance.

The complete TheraBite ActiveBand Kit includes the following items (one each):

ActiveBand (Fig. 1a)

Load Scale (Fig. 1b)

Instructions for use

Exercise Log

1.4 Expectations of the Device

The TheraBite ActiveBand is an accessory for active training of the jaw closing muscles. It can be used as a stand-alone treatment method, or as a complement to passive motion training.

It can be used to improve muscle strength and endurance in weakened muscles, and it can also be used to prevent decreased muscle strength during (chemo) radiation in head and neck cancer patients. The level of resistance can be adjusted from approximately 5N to 60N, (5 N equals about 1/2 kilogram of resistance).

NOTE: The Load Scale is only an indication of the approximate resistance measured in Newton (N), and is intended for reference only.

1.5 WARNINGS

DO NOT

- Do not use the ActiveBand before consulting a clinician if you have a specific neuro-muscular diseases, for example ALS or similar.
- Do not continue to exercise if you feel severe pain at any time while exercising with the ActiveBand. You should immediately stop using the product and contact your clinician.

- Do not apply excessive force or overexert yourself, this may cause injury to the jaw muscles and the temporomandibular joint* (TMJ) region.

*The joint where the lower jaw (the mandible) joins the temporal bone of the skull. It is located immediately in front of the ear on each side of the head.

2. Instructions for Use

2.1 Preparation

Before each use check that the device is undamaged. Do not use the device if it is damaged.

Attaching the Load Scale (Fig. 2a-c)

1. Make sure the surface of the Jaw Mobilizer is clean and dry before attaching the Load Scale.
2. Peel off the liner from the Load Scale (Fig. 2a).
3. Attach the Load Scale to the bottom of the Jaw Mobilizer as shown in (Fig. 2b and 2c).

NOTE: Load resistance 1 should point towards the Adjustment Knob and Load resistance 6 towards the Mouthpieces.

Attaching the ActiveBand (Fig. 3a-f)

CAUTION: Do not use lubricating agents on the TheraBite ActiveBand. The ActiveBand may change position during training and/or change its properties causing it to break.

1. Ensure that the Fine Adjustment Knob on the TheraBite is screwed in as far as possible.
2. Loop the ActiveBand around the Jaw Mobilizer from the tip of the device close to the Adjustment Knob (Fig. 3a and 3b).
3. To move the ActiveBand on the Load Scale, loosen the knob until the Mouthpieces are at maximum opening. Slide the ActiveBand along the Load Scale in order to align it to one of the bars that are illustrated beside every load number (1–6), (Fig. 3c). Number 1 means very low resistance and number 6 is very high. Fig. 3d illustrates ActiveBand aligned at position 3. Make sure that the ActiveBand is in a vertical position.
4. Adjust the Fine Adjustment Knob while applying pressure to the Lever so that the lower edge of the Mouthpiece opens to approximately 15 mm, or any opening recommended by your clinician (Fig. 3e and 3f).

NOTE: An opening greater than 25 mm is not recommended.

CAUTION: Do not overstretch or leave the ActiveBand in stretched position when not in use, this may affect the durability and strength of the ActiveBand.

2.2 Rehabilitation program (Fig. 4a-c)

It is recommended that your clinician selects and combines exercises that are appropriate for you, including regular monitoring and review.

1. Hold the Jaw Mobilizer in your hand.
2. Press the Mouthpieces together with your fingers to insert, and bite very gently (Fig. 4a).
3. Hold position, open and then close your mouth against the resistance of the ActiveBand (Fig. 4b and 4c).

CAUTION: Do not release the Lever abruptly as the Mouthpiece is preloaded with the force applied by the ActiveBand. Release the Lever slowly to avoid pain or injury to your teeth or jaw.

4. Hold for one second, and then repeat the exercise (Fig. 4b and 4c).

The position of the ActiveBand on the Load Scale, number of repetitions, number of treatment sessions per day, and length of treatment time will vary depending on your individual situation, goals and current muscle strength. Your clinician will advise a schedule based on your individual situation.

For example, to increase muscle strength, a 3-5-30 schedule where you perform the exercises 3 times daily, using 5 sets of 30 repetitions separated by one minute rest between each set can be used.

To maintain muscle strength (for example during (chemo) radiotherapy), a less intensive 2-4-10 schedule where you perform the exercise 2 times per day, doing 4 sets of 10 repetitions can be used.

In many cases, these exercises are combined with passive stretching exercises as described in the TheraBite instructions for use.

NOTE: Mild fatigue may be experienced in the early stages of exercise, but you should not experience excessive fatigue. Excessive fatigue may be due exercising too long or with a too high resistance. Take a break from exercising if you experience excessive fatigue. If necessary, contact your clinician.

Too long exercise at too high resistance or too large mouth opening can also lead to pain in the joint. If this occurs reduce the quantity of exercises, resistance and/or mouth opening. If necessary, contact your clinician.

Start at a resistance where you can easily do 30 repetitions, and then gradually increase the resistance, following the guidance of your clinician.

Your minimum exercise plan: The resistance at which you can comfortably carry out your prescribed exercise routine.

Your maximum exercise plan: The maximum resistance at which you can perform the same sets of exercises.

2.3 Use of Exercise Log

Use the Exercise Log to document and monitor your progress. The columns in the log are for the following data:

- A: Date
- B: Resistance
- C: Frequency of exercise / day
- D: Number of sets
- E: Number of repetitions

The resistance (B above) on the load scale goes from 1 to 6.

Scale	Resistance
1	Very low
2	Low
3	Medium low
4	Medium high
5	High
6	Very high

2.4 Cleaning

Use only water and dish soap to clean the ActiveBand. Rinse well, shake off excess water, and dry.

2.5 Device lifetime

Six months of daily use.

2.6 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Coloplast A/S. All rights reserved.
Atos and the Atos Medical logo are trademarks of Coloplast A/S.
Atos Medical UK Ltd Tottle Road
Atos Medical AB